



Ultrasonic gel



All professional ultrasound scanners will use ultrasound gel at time of scanning. The primary purpose of this acoustic couplant is to allow the ultrasound waves to pass into the body without reflecting off intervening air. To further optimise their images, many scanners will also ask to shave the bitch underneath. This enables even better contact between the probe and the skin, as even a light covering of hair will trap air

ADVANTAGES

- 1.No dry skin
- 2.High durability
- 3.Without bubbles
- 4.Without fatty and alcohol
- 5.With high conductivity
- 6.Without paraben and bad smell
- 7.In volumes of 260 cc and 5 liters

